

NORMA Face Shield

Information from manufacturer NORMA Distribution Center GmbH
As per PPE Regulation (EU) 2016/425
Version: 28.05.2020

Manufacturer: NORMA Distribution Center GmbH, Unterm Ohmberg 24, 34431 Marsberg, Germany



Instructions for

1. Storage: Store the NORMA Face Shield at room temperature in a dry, closed space.

2. Use: Before placing the NORMA visor on the NORMA frame, make sure to remove the protective film, which is applied to both sides of the shield. The protective film is coloured green on one side and is clear on the other. Tip: A piece of tape helps to remove the protective film. Simply apply it to the edge of the film so that the piece of tape sticks out a bit. Then, pull it off. The visor is clipped to the band by placing the holes onto the pins. The top rear ends of the visor are inserted into the small side bands. The visor does not have to be adjusted to the user.

We recommend using the product described above once. However, if cleaned and disinfected properly, the product can also be used multiple times. As a precaution, we would like to point out that this recommendation is non-binding, non-validated and non-certified. The user is responsible for verifying which processes and measures are suitable and necessary for cleaning and disinfecting the product for repeat use in order to meet or achieve the properties required for the intended use of the product. This applies in particular when our products are processed, combined or used otherwise.

3. Cleaning and disinfection: The visor can be cleaned using a microfibre cloth. It can be disinfected using customary agents.

4. Maintenance and inspection: Inspect the PPE before use. We recommend replacing defective PPE (e.g. scratched visor).

PPE functions:

1. This eye protector protects against liquid droplets and splashes, reducing the risk of a droplet infection. Due to its design, it does not provide full protection against viruses, bacteria and other biological agents.
2. The product does not provide protection against mechanical actions.
3. The product can be combined with face masks and safety eyewear.
4. Dispose of the product after use! Re-use is up to the user's discretion.

Accessories:

Frame and visor are available separately from your sales partner. The relevant NORMA parts numbers are:
NORMA Frame: 0700 8701 000; NORMA Visor: 0700 8702 000.

PPE expiration date:

The PPE has a maximum storage life of 5 years.

Packaging type:

No special packaging type required for transport;
the visor is protected by protective film on both sides.

Meaning of the labels:

NDC CE

1. NDC = Manufacturer abbreviation for
NORMA Distribution Center GmbH
2. CE = CE mark

Notified body:

DIN CERTCO Gesellschaft für Konformitätsbewertung
mbH, Alboinstraße 56, 12103 Berlin, Germany,
(Notified Body No.: 0196).

Harmonisation legislation:

1. Regulation (EU) 2016/425
2. Directive 2001/95/EC
3. EN 166:2001

Website for EU declaration of conformity:

http://www.normagroup.com/norma.nsf/id/Landing-NORMA-Face-Shield_EN

All information from the manufacturer risk assessment:

1. This eye protector does not provide protection against mechanical actions.
2. This eye protector does not provide protection against radiation (e.g. solar radiation, UV, welding, IR).
3. This eye protector may only be used at room temperature.
4. This eye protector is not fire-resistant.
5. Sharp edges can result in injuries.
6. Do not use a damaged eye protector, as this can impair your vision and protection.
7. If the eye protector is damaged, it must be replaced.
8. Dispose of the eye protector after use. Re-use is up to the user's discretion.
9. This eye protector does not provide protection against aggressive fluids.

NORMA Distribution Center GmbH



How do you like our NORMA Face Shield?
Do you have any questions or comments?

We are happy to receive your feedback: norma.to/faceshield

